



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Mouse anti-Human CD62P, clone AK-6 (Monoclonal)

Clone no. AK-6

MONOSAN

---

Product name	Mouse anti-Human CD62P, clone AK-6 (Monoclonal)
Host	Mouse
Applications	IHC-fr, FC, IF, IHC-P, WB
Species reactivity	human
Conjugate	-
Immunogen	Human platelet membrane glycoproteins
Isotype	IgG1
Clonality	Monoclonal
Clone number	AK-6
Size	0.2 mg
Concentration	IgG 1 mg/ml
Format	Purified
Storage buffer	PBS with <0.1% sodium azide
Storage until expiry date	aliquots -20°C. Thawed 2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Human CD62P, clone AK-6 (Monoclonal)

Clone no. AK-6

MONOSAN

**Additional info**

Mouse anti Human CD62P antibody, clone AK-6 recognizes the CD62P, also known as P-selectin, Granule membrane protein 140, GMP140, Leukocyte-endothelial cell adhesion molecule 3 or Platelet activation dependent granule-external membrane protein. CD62P is a 830 amino acid, including a 41 amino acid signal peptide, ~140 kDa single pass type I transmembrane glycoprotein expressed on activated platelets and endothelial cell

**References**

1. Dunlop LC et al. J Exp Med 1992; 175: 1147-50
2. Sopper S et al. Cytometry 1997; 29: 351-62
3. Knipe L et al. Blood 2010; 116: 2183-91
4. Turner CP et al. Br J Haematol 2003; 121: 623-31
5. Cipok M et al. Eur J Haematol 2019; 103: 225-233

**FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES**