



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Mouse anti-H Type 2, clone 19-OLE (Monoclonal)

Clone no. 19-OLE

MONOSAN

Product name	Mouse anti-H Type 2, clone 19-OLE (Monoclonal)
Host	Mouse
Applications	FC, IHC-fr, IHC-P, IF, ICC
Species reactivity	human
Conjugate	-
Immunogen	mucin isolated from a mucinous colonic adenocarcinoma
Isotype	IgM
Clonality	Monoclonal
Clone number	19-OLE
Size	100 ug
Concentration	100 ug/ml
Format	-
Storage buffer	PBS with 0.02% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-H Type 2, clone 19-OLE (Monoclonal)

Clone no. 19-OLE

MONOSAN

Additional info

19-OLE reacts with H type 2 antigen, the basis of the ABO blood group system, involving three carbohydrate antigens: A, B, and H. A, B, and AB individuals express a glycosyltransferase activity that converts the H antigen to the A antigen (by addition of UDP-GalNAc) or to the B antigen (by addition of UDP-Gal), whereas O individuals lack such activity. It is expressed on endothelial cells, epithelial cells and granulocytes. Increased expression of this antigen has been observed on some tumor tissues such as gastric carcinomas, urothelial carcinomas, and colon carcinomas.

References

1. Bara, J. et al. Blood transfusion and immunohaematology, Ph Rouger, D Anstee
2. Marionneau, S. et al. Gasterenterology, 122: 1967-1877 (2002)
3. Garcher, K. et al. Invest Ophthalmol Vis Sci. 35: 1184-1191 (1994)
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES