



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Product datasheet      MON3206

**MONOSAN**<sup>®</sup>

Mouse anti-bcl-2, clone 124 (Monoclonal)

Clone no.      124

MONOSAN

---

Product name	Mouse anti-bcl-2, clone 124 (Monoclonal)
Host	Mouse
Applications	IHC-P (1:100-1:500)
Species reactivity	human
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	IgG1-k
Clonality	Monoclonal
Clone number	124
Size	1 ml
Concentration	n/a
Format	-
Storage buffer	Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-bcl-2, clone 124 (Monoclonal)

Clone no. 124

MONOSAN

---

**Additional info**

Bcl-2 is the best characterized protein family involved in regulation of apoptotic cell death, consisting of anti-apoptotic and pro-apoptotic members. Bcl-2 is a useful marker for identifying neoplastic cells in follicular lymphoma. Antibodies specific for the Bcl-2 protein can be used to distinguish between reactive and neoplastic follicular proliferation in lymph node biopsies.

**References**

1. Tsujimoto, Y. Genes Cells. 1998; 3:697-707
2. Gaulard P, et al. Am J Pathol. 1992; 140: 1089-95
3. Wang T, et al. APMIS. 1995; 103:655-62
4. West RB, et al. Am J Clin Pathol. 2002; 117:636-43
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES