



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Mouse anti-CD25, clone 4C9 (Monoclonal)

Clone no. 4C9

MONOSAN

| | |
|---------------------------|---|
| Product name | Mouse anti-CD25, clone 4C9 (Monoclonal) |
| Host | Mouse |
| Applications | IHC-P (1:10-1:50) |
| Species reactivity | human |
| Conjugate | - |
| Immunogen | Unknown or proprietary to MONOSAN and/or its suppliers |
| Isotype | IgG2b |
| Clonality | Monoclonal |
| Clone number | 4C9 |
| Size | 1 ml |
| Concentration | n/a |
| Format | - |
| Storage buffer | Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide |
| Storage until expiry date | 2-8°C |

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-CD25, clone 4C9 (Monoclonal)

Clone no. 4C9

MONOSAN

Additional info

CD25, Interleukin-2 receptor alpha chain, is the alpha subunit of the cell surface receptor which regulates regulatory T-cells. CD25 has been detected in various hematological malignancies including adult T-cell leukemia/lymphoma, and hairy cell leukemia. MON 3230 has also been useful in identifying mast cells in skin biopsies in the setting of Urticaria Pigmentosa, which is predictive of Systemic Mastocytosis.

References

1. Hahn HP, et al. Am J Surg Pathol. 2007 Nov;31(11):1669-76
2. Hollmann TJ, et al. Am J Surg Pathol. 2008 Jan;32(1):139-45
3. L  tourneau S, et al. Clin Immunol. 2009; 123:758-62
4. De Toter D, et al.. Leuk Lymphoma. 1994; 104:412-9
5. Qayyum S, et al. Archives of Pathology & Lab Med. 2014; 138:282-6

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES