



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Mouse anti-CEA, clone CEA31 (Monoclonal)

Clone no. CEA31

MONOSAN

| | |
|---------------------------|---|
| Product name | Mouse anti-CEA, clone CEA31 (Monoclonal) |
| Host | Mouse |
| Applications | IHC-P (1:100-1:500) |
| Species reactivity | human |
| Conjugate | - |
| Immunogen | Unknown or proprietary to MONOSAN and/or its suppliers |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| Clone number | CEA31 |
| Size | 1 ml |
| Concentration | n/a |
| Format | - |
| Storage buffer | Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide |
| Storage until expiry date | 2-8°C |

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-CEA, clone CEA31 (Monoclonal)

Clone no. CEA31

MONOSAN

Additional info

Anti-CEA is employed as a tool to assist in the distinction between adenocarcinoma and mesotheliomas, along with other markers such as calretinin, CK 5/6, D2-40, HBME-1, Napsin A, MOC31, and Ber-EP4. Anti-CEA positivity is seen in adenocarcinomas from the lung and colon, as well.

References

1. Go, VLW, et al., Cancer 1976;37:562-566
2. Delellis, RA, et al., Am J Clin Pathol 1978;50:587-594
3. Abutaily AS et al. J Clin Pathol. 2002 Sep;55(9):662-8
4. Bhatnagar J et al. Anticancer Res. 2002;22(3):1849-57
5. Carella R. et al. Am J Surg Pathol. 2001 Jan;25(1):43-50

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES