



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Rabbit anti-Human GFAP, clone EP672Y (Monoclonal)

Clone no. EP672Y

MONOSAN

| | |
|---------------------------|---|
| Product name | Rabbit anti-Human GFAP, clone EP672Y (Monoclonal) |
| Host | Rabbit |
| Applications | IHC-P (1:25-1:100) |
| Species reactivity | human |
| Conjugate | - |
| Immunogen | Unknown or proprietary to MONOSAN and/or its suppliers |
| Isotype | IgG |
| Clonality | Monoclonal |
| Clone number | EP672Y |
| Size | 1 ml |
| Concentration | n/a |
| Format | - |
| Storage buffer | Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide |
| Storage until expiry date | 2-8°C |

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Rabbit anti-Human GFAP, clone EP672Y (Monoclonal)

Clone no. EP672Y

MONOSAN

Additional info

The antibody detects astrocytes, Schwann cells, satellite cells, enteric glial cells, and some groups of ependymal cells. This marker is mainly used to distinguish neoplasms of astrocytic origin from other neoplasms in the central nervous system.

References

1. Jessen, KR, et al. J Neurosci 1983;3:2206-2218
2. Regner A et al. J Neurotrauma. 2001 Aug;18(8):783-92
3. Nagashima G et al. Clin Neurol Neurosurg. 2002 May;104(2):125-31
4. Choi, BH, et al. Science 1984;223:407-409
5. Viale, G, et al. Virchow's Arch A Pathol Anat 1991;418:339-348

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES