

Produktinformation



Forschungsprodukte & Biochemikalien
Zellkultur & Verbrauchsmaterial
Diagnostik & molekulare Diagnostik
Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



Lieferung & Zahlungsart siehe unsere Liefer- und Versandbedingungen

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien T. +43(0)1 489 3961-0 F. +43(0)1 489 3961-7 <u>mail@szabo-scandic.com</u> www.szabo-scandic.com

Product datasheet

MON4009



Mouse anti-Alpha-Tubulin, clone TU-01 (Monoclonal) Clone no. TU-01

| Product name | Mouse anti-Alpha-Tubulin, clone TU-01 (Monoclonal) |
|---------------------------|-----------------------------------------------------------------------------------------------------|
| Host | Mouse |
| Applications | WB, IHC-P, ICC, IP, FC |
| Species reactivity | Mouse, Pig, Human, Turkey, Eisenia, Paramecium, Nicotiana, Yeast, Arabidopsi: |
| Conjugate | - |
| Immunogen | Fraction of tubulin purified from porcine brain by two cycles of polymerization - depolymerization. |
| lsotype | lgG1 |
| Clonality | Monoclonal |
| Clone number | TU-01 |
| Size | 0.1 mg |
| Concentration | 1 mg/ml |
| Format | - |
| Storage buffer | Phosphate buffered saline (PBS) solution with 15 mM sodium azide |
| Storage until expiry date | 2-8°C |

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Product datasheet

MON4009

MONOSAN

Mouse anti-Alpha-Tubulin, clone TU-01 (Monoclonal) Clone no. TU-01

MONOSAN

Additional info The antibody TU-01 recognizes a defined epitope (aa 65-97) on N-terminal structural domain of alpha-tubulin. The microtubules are intracellular dynamic polymers made up of evolutionarily conserved polymorphic alpha/beta-tubulin heterodimers and a large number of microtubuleassociated proteins (MAPs). The microtubules consist of 13 protofilaments and have an outer diameter 25 nm. Microtubules have their intrinsic polarity; highly dynamic plus ends and less dynamic minus ends. Microtubules are required for vital processes in eukaryotic cells including mitosis, meiosis, maintenance of cell shape and intracellular transport. Microtubules are also necessary for movement of cells by means of flagella and cilia. In mammalian tissue culture cells microtubules have their minus ends anchored in microtubule organizing centers (MTOCs). The GTP (guanosintriphosphate) molecule is an essential for tubulin heterodimer to associate with other heterodimers to form microtubule. In vivo, microtubule dynamics vary considerably. Microtubule polymerization is reversible and a populations of microtubules in cells are on their minus ends either growing or shortening – this phenomenon is called dynamic instability of microtubules. On a practical level, microtubules can easily be stabilized by the addition of nonhydrolysable analogues of GTP (eg. GMPPCP) or more commonly by anticancer drugs such as Taxol. Taxol stabilizes microtubules at room temperature for many hours. Using limited proteolysis by enzymes both tubulin subunits can be divided into N-terminal and C-terminal structural domains. The alpha-tubulin (relative molecular weight around 50 kDa) is globular protein that exists in cells as part of soluble alpha/beta-tubulin dimer or it is polymerized into microtubules. In different species it is coded by multiple tubulin genes that form tubulin classes (in human 6 genes). Expressed tubulin genes are named tubulin isotypes. Some of the tubulin isotypes are expressed ubiquitously, while some have more restricted tissue expression. Alpha-tubulin is also subject of numerous post-translational modifications. Tubulin isotypes and their posttranslational modifications are responsible for multiple tubulin charge variants - tubulin isoforms. References 1. Viklicky V, et al. Cell Biol Int Rep. 1982 Aug;6(8):725-31 2 Grimm M, et al. Biochim Biophys Acta. 1987 Jul 24;914(1):83-8. 3. Linhartova I, et al. Biochem J. 1992 Dec 15;288 (Pt 3):919-24. 4. Draber P, et al. Eur J Cell Biol. 1986 Jun;41(1):82-8 5.

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

www.monosan.com