



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Mouse anti-E-Selectin CD62-E, clone ENA2, Purified (Monoclonal)

Clone no. ENA2

MONOSAN

Product name	Mouse anti-E-Selectin CD62-E, clone ENA2, Purified (Monoclonal)
Host	Mouse
Applications	IHC-fr,FUNC,ELISA,IF
Species reactivity	human, baboon, chimpanzee, cynomolgus monkey, rhesus monkey
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	IgG1
Clonality	Monoclonal
Clone number	ENA2
Size	1 ml
Concentration	100 ug/ ml
Format	-
Storage buffer	PBS with 0.1% BSA and 0.02% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-E-Selectin CD62-E, clone ENA2, Purified (Monoclonal)

Clone no. ENA2

MONOSAN

Additional info

ENA2 reacts with E-selectin CD62-E, previous designated the Endothelial Leucocyte Adhesion Molecule-1 (ELAM-1). The antibody reacts with human endothelial cells activated with TNF-alpha, IL-1 or endotoxin. The antibody was found to react also with cells transfected with the E-selectin gene. The antibody inhibits the adhesion of granulocytes both neutrophilic and eosinophilic.

References

1. Leeuwenberg; JFM et al. Eur J Immunol 1989; 19: 715
2. Leeuwenberg, JFM et al Clin and Exp Immunol 1990, 81: 496
3. Leeuwenberg; JFM et al. Immunology 1990; 71: 301
4. Leeuwenberg JFM et al. J Immunology 1990; 145: 2110
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES