



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Mouse anti-Hairy Cell Leukemia (HCL), clone DBA.44 (Monoclonal)

Clone no. DBA.44

MONOSAN

---

|                           |   |
|---------------------------|---|
| Product name              | Mouse anti-Hairy Cell Leukemia (HCL), clone DBA.44 (Monoclonal) |
| Host                      | Mouse   |
| Applications              | IHC-P (1:15-1:30)   |
| Species reactivity        | human   |
| Conjugate                 | -   |
| Immunogen                 | DEAU cell line established from a large cell lymphoma           |
| Isotype                   | IgMk  |
| Clonality                 | Monoclonal  |
| Clone number              | DBA.44  |
| Size                      | 1 ml  |
| Concentration             | lot specific  |
| Format                    | -   |
| Storage buffer            | Tissue culture supernatant with sodium azide                    |
| Storage until expiry date | 2-8°C   |

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Hairy Cell Leukemia (HCL), clone DBA.44 (Monoclonal)

Clone no. DBA.44

MONOSAN

**Additional info**

This antibody reacts with a subset of B-lymphocytes localized in the follicular mantle zone. It reacts with 97% of hairy cell leukemia cases. This antibody shows strong positive staining of about 35% of cases of high grade B cell lymphomas. Positive control Tonsil.

**References**

1. Al Saati et al. 1989;74:2476
2. -
3. -
4. -
5. -

**FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES**