



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Mouse anti-Melanoma Associated Antigen, clone NKI/C3 (Monoclonal)

Clone no. NKI/C3

MONOSAN

| | |
|---------------------------|---|
| Product name | Mouse anti-Melanoma Associated Antigen, clone NKI/C3 (Monoclonal) |
| Host | Mouse |
| Applications | IHC-P (1:10-1:20) |
| Species reactivity | human |
| Conjugate | - |
| Immunogen | Unknown or proprietary to MONOSAN and/or its suppliers |
| Isotype | IgG1 |
| Clonality | Monoclonal |
| Clone number | NKI/C3 |
| Size | 1 ml |
| Concentration | n/a |
| Format | - |
| Storage buffer | 50 mM PBS at pH 7.4 with 1% BSA and 15mM sodium azide |
| Storage until expiry date | 2-8°C |

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Melanoma Associated Antigen, clone NKI/C3 (Monoclonal)

Clone no. NKI/C3

MONOSAN

Additional info

This antibody recognizes a heterogeneous 25-110 kD glycoprotein that is located mainly in the inner side of membranes of cytoplasmic vesicles in melanoma cells. The antigen has a 25 kD unglycosylated precursor in formalin fixed and paraffin-embedded tissue sections. Cross reactivity: The antibody reacts with some mucus producing tumors, carcinoids, carcinomas of the thyroid, mast cells, histiocytes in tumor regions and with cells with secretory functions such as salivary glands, bronchial glands, sweat glands, pancreas and prostate. The antibody should be used on formalin fixed paraffin embedded tissue sections, it is not advisable to use the antibody on frozen sections.

References

1. MacKie, R.M., et al., J. Clin. Pathol. 37, 367 (1984)
2. van Duinen, S.G., et al., Cancer 53, 1566 (1984)
3. Vennegoor, C., et al., Int. J. Cancer 35, 287 (1985)
4. Palmer, A.A., et al., Pathology 17, 335 (1985)
5. Hagen, E.C., et al., Histopathology 10, 689 (1986)

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES