



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

## Mouse anti-Herpes simplex virus type 1, clone T96 (Monoclonal)

Clone no. T96

MONOSAN

---

Product name	Mouse anti-Herpes simplex virus type 1, clone T96 (Monoclonal)
Host	Mouse
Applications	WB, ICC, ELISA
Species reactivity	Herpes Simplex Virus
Conjugate	-
Immunogen	Extract of HSV-1 infected VERO (green monkey kidney) cells.
Isotype	IgG2b
Clonality	Monoclonal
Clone number	T96
Size	0.1 mg
Concentration	1 mg/ml
Format	-
Storage buffer	Phosphate buffered saline (PBS) solution with 15 mM sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

## Mouse anti-Herpes simplex virus type 1, clone T96 (Monoclonal)

Clone no. T96

MONOSAN

**Additional info**

Herpes simplex Vvrus (HSV) is a virus that manifests itself in two common viral infections. There are actually two types of herpes simplex virus, HSV1 and HSV2. These are very similar in many ways, and both can cause either oral herpes or genital herpes. HSV1 - most commonly develops into oral herpes infecting the lips (fever blisters or cold sores). HSV1 can also infect the genital area causing sores to develop. HSV2 - generally infects the genital area (genital herpes); however, HSV2 can also infect the mouth.

**References**

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES