



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Mouse anti-Leukemia cell line HL-60, clone IPO-M6 (Monoclonal)

Clone no. IPO-M6

MONOSAN

Product name	Mouse anti-Leukemia cell line HL-60, clone IPO-M6 (Monoclonal)
Host	Mouse
Applications	FC, IHC-fr, IP
Species reactivity	human
Conjugate	-
Immunogen	human cell line HL-60
Isotype	IgG1
Clonality	Monoclonal
Clone number	IPO-M6
Size	100 ug
Concentration	100 ug/ml
Format	-
Storage buffer	PBS with 0.02% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Leukemia cell line HL-60, clone IPO-M6 (Monoclonal)

Clone no. IPO-M6

MONOSAN

Additional info

IPO-M6 reacts with human leukemia cell line HL-60 and immuno-precipitates two proteins with MW of 48 and 52 kDa. IPO-M6 does not stain B cell lines Daudy, PHS, Namalwa, RPMI-1788 and T-cell lines CCRF-HSB2, Jurkat and Molt-4. IPO-M6 can be applied for staining of monocytes and up to 10 % of lymphocytes from peripheral blood of healthy donors. Blast cells of patients with AMMonL (M5 following FAB classification), AMMonL (M4) and hairy cells leukemia are IPO-M6 positive. The antigen, defined by IPO-M6, is particularly expressed on blood cells from patients with infectious mononucleosis and CLL. Histiocytes and macrophages are also positive. Malignant cells from patients with AML (M1 and M2), T-ALL, B ALL are IPO-M6 negative.

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES