



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Mouse anti-HLA-B27, clone EP-4 (Monoclonal)

Clone no. EP-4

MONOSAN

Product name	Mouse anti-HLA-B27, clone EP-4 (Monoclonal)
Host	Mouse
Applications	IHC-fr, FC, IF
Species reactivity	human
Conjugate	-
Immunogen	Lymphocytes from HLA-B27 individual
Isotype	IgM,kappa
Clonality	Monoclonal
Clone number	EP-4
Size	100 ug
Concentration	100 ug/ml
Format	-
Storage buffer	PBS with 0.02% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-HLA-B27, clone EP-4 (Monoclonal)

Clone no. EP-4

MONOSAN

Additional info

EP-4 recognizes the HLA-B27 cell surface antigen on human cells. HLA-B27 has been found to be highly associated (60%) with ankylosing spondylitis (Bekhterev's disease, MarieStrümpell disease or "bamboo spine"). While rheumatoid factor tests will be negative, testing for the presence of HLA-B27 in a patient can confirm the diagnosis of ankylosing spondylitis. Also postgonococcal arthritis and acute anterior uveitis are associated with HLA-B27.

References

1. Hansen JA et al. Hematol Oncol Clin North Am, 4(3): 507-515 (1990)
2. El-Shabrawi, Y. et al. Ophthalmology 113: 695-700 (2006).
3. Kellner H, Clin Investig. 72(4): 321-8 (1994)
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES