



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Mouse anti-MRP1, clone MRPm6 (Monoclonal)

Clone no. MRPm6

MONOSAN

Product name	Mouse anti-MRP1, clone MRPm6 (Monoclonal)
Host	Mouse
Applications	FC, WB, ICC, IHC-fr, IHC-P (1:20-1:50)
Species reactivity	human
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	IgG1
Clonality	Monoclonal
Clone number	MRPm6
Size	1 ml
Concentration	250 ul/ ml
Format	-
Storage buffer	Serum free tissue culture supernatant with 0.7% BSA and 0.1% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-MRP1, clone MRPm6 (Monoclonal)

Clone no. MRPm6

MONOSAN

Additional info

MRPm6 reacts with an internal epitope of MRP, a 180-195 kD transmembrane transporter protein overexpressed in various human non-P-glycoprotein MDR tumor cell lines. MRPm6 was raised against a bacterial fusion protein of MRP, containing a segment of 170 amino acids in the carboxy terminal end and part of the carboxy proximal nucleotide binding domain of the protein. MRPm6 does not crossreact with the human MDR1 and MDR3 gene products

References

1. Moll I et al. Eur J Cell Biol 2005; 84: 259-271
2. Riedel I et al. Virchows Arch 2001; 438: 181-191
3. Romih R et al. Cell Biol 1998; 109: 263-269
4. Demirkesen C et al. J Cutan Pathol. 1995; 22: 518-535
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES