



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Product datasheet      MON9200

**MONOSAN**<sup>®</sup>

MDR Sampler Pack

Clone no.      Sampler

MONOSAN

---

Product name	MDR Sampler Pack
Host	Rat
Applications	ICC, IHC-fr (1:20), IHC -P, WB, FC
Species reactivity	human
Conjugate	-
Immunogen	human MRP1 amino-proximal half segment of 168 aa
Isotype	IgG2a
Clonality	Monoclonal
Clone number	Sampler
Size	250ul
Concentration	250 ul/ ml
Format	-
Storage buffer	cell culture supernatant with 0.7% BSA and 0.1% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

---

**Additional info**

MRPr1 reacts with an epitope of MRP, a 180-195 kD transmembrane transporter protein overexpressed in various human non-P-glycoprotein MDR tumor cell lines. MRPr1 was raised against a bacterial fusion protein of MRP, containing a segment of 168 amino acids in the amino-proximal half of the protein. MRPr1 does not cross-react with the human MDR1 and MDR3 gene products (Flens et al. 1994).

**References**

1. Flens M et al. Cancer Res 1994; 54: 4557-4564
2. Cole S et al. Science 1992; 258: 1650-1654
3. Zaman et al. Proc Natl Acad Sci 1994; 91: 8822-8826
4. Flens M et al. Am J Pathol 1996
5. Schandoreff et al Am J Pathol 1995

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES