



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Mouse anti-Blood Group A (CD173), clone HE-10 (monoclonal)

Clone no. HE-10

MONOSAN

Product name	Mouse anti-Blood Group A (CD173), clone HE-10 (monoclonal)
Host	Mouse
Applications	FC, IHC-fr, IHC-P, IF, FUNC
Species reactivity	human
Conjugate	-
Immunogen	mixture of erythrocytes of blood group A1 and a glycoprotein fraction isolated from the saliva of secretors with blood group A of human origin
Isotype	IgM-K
Clonality	Monoclonal
Clone number	HE-10
Size	100 ug
Concentration	100 ug/ml
Format	-
Storage buffer	PBS with 0.02% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Blood Group A (CD173), clone HE-10 (monoclonal)

Clone no. HE-10

MONOSAN

Additional info

HE-10 preferably reacts with determinants of chain A type 3 and 4 and chain H type 3 and 4, but not with type 1 and 2 chain structures. It is not reactive with immunodominant A trisaccharide. Increased expression of this antigen has been observed on some tumor tissues such as gastric carcinomas, urothelial carcinomas, and colon carcinomas. HE-10 does not react with normal tissue sections of donors with blood group B and 0 but it reacts specifically with malignant tissues in these individuals. HE-10 is applicable for red cell agglutination, tissue staining and immunofluorescence tests.

References

1. Němec M. et al. Vox Sang 52:125-8 (1987)
2. Vanák J. et al. Neoplasma 36(4): 479-487 (1989)
3. Tichý, M. et al. Neoplasma 37(4): 451-459 (1990)
4. Tichý, M. et al. Acta Univ Palacki Olomuc Fac Med. 126: 57-69 (1990)
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES