



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Product datasheet MON9912

MONOSAN[®]

Mouse anti-MALT1, clone EBS-T-007 (Monoclonal)

Clone no. EBS-T-007

MONOSAN

Product name Mouse anti-MALT1, clone EBS-T-007 (Monoclonal)

Host Mouse

Applications IHC-fr, IHC-P, WB

Species reactivity human

Conjugate -

Immunogen human MALT1 recombinant protein

Isotype IgG1-K

Clonality Monoclonal

Clone number EBS-T-007

Size 100 ug

Concentration 100 ug/ml

Format -

Storage buffer PBS with 0.02% sodium azide

Storage until expiry date 2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-MALT1, clone EBS-T-007 (Monoclonal)

Clone no. EBS-T-007

MONOSAN

Additional info

EBS-T-007 reacts with paracaspase MALT1, crucial for B- and T-cell activation/proliferation and activation of transcription factor NF- κ B. It has an N-terminal death domain, two central immunoglobulin-like domains involved in the binding to the B-cell lymphoma 10 (BCL10) protein and a caspase-like domain. MALT1 and BCL10, can both be found translocated in MALT lymphoma. In such cases either BCL10 or MALT1 or both are highly expressed, depending on the site of translocation. Normal cells and MALT lymphomas lacking translocations exhibit much lower levels of expression.

References

1. Ye H, et al, Am J Pathol 157: 1147-54 (2000)
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES