



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Mouse anti-CD55, clone F4-29D9 (Monoclonal)

Clone no. F4-29D9

MONOSAN

Product name	Mouse anti-CD55, clone F4-29D9 (Monoclonal)
Host	Mouse
Applications	FC, IHC-fr, IF
Species reactivity	human
Conjugate	-
Immunogen	human umbilical vein endothelial cells
Isotype	IgG1-K
Clonality	Monoclonal
Clone number	F4-29D9
Size	100 ug
Concentration	100 ug/ml
Format	-
Storage buffer	PBS with 0.02% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-CD55, clone F4-29D9 (Monoclonal)

Clone no. F4-29D9

MONOSAN

Additional info

F4-29D9 reacts with CD55 or DAF (Decay Accelerating Factor). All leucocytes as well as human erythrocytes, fibroblasts, platelets, endothelial cells and neuroectodermal cells are positive for DAF. F4-29D9 also recognizes an antigen on spermatozoid cells. It is a glycosylphosphatidylinositol anchored (GPI-anchored) member of the membrane bound complement regulatory proteins that inhibit autologous complement cascade activation. CD55 also serves as receptor for CD97 and for echovirus and coxsackie B virus. DAF is deficient in both granulocytes and monocytes in patients with paroxysmal nocturnal haemoglobinuria.

References

1. Kishimoto T. et al., eds. Leukocyte Typing VI, p509-514, Garland Publishing, Inc.
2. B.E. Loveland - in Leukocyte Typing VI - Part 6 CD55 Workshop Panel Report p1
3. Ruix-Delgado, GJ et al, Hematology 14: 33-7 (2009)
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES