



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Rat anti-Mouse CD8a (mouse), clone IBL-3/25 (Monoclonal)

Clone no. IBL-3/25

MONOSAN

---

Product name	Rat anti-Mouse CD8a (mouse), clone IBL-3/25 (Monoclonal)
Host	Mouse
Applications	FC, IHC-fr, IF
Species reactivity	mouse
Conjugate	-
Immunogen	mouse lymph node cells
Isotype	IgG-K
Clonality	Monoclonal
Clone number	IBL-3/25
Size	100 ug
Concentration	100 ug/ml
Format	-
Storage buffer	PBS with 0.02% sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Rat anti-Mouse CD8a (mouse), clone IBL-3/25 (Monoclonal)

Clone no. IBL-3/25

MONOSAN

---

**Additional info**

IBL-3/25 is directed against the  $\alpha$ -chain of CD8. The CD8 complex consists of a disulphide-linked  $\alpha/\beta$  heterodimer with a MW of 30 kDa or an  $\alpha/\alpha$  homodimer with a MW of 32 kDa. The CD8 molecule binds to HLA class I molecules during interaction of CD8+ T-cells with antigenpresenting cells or with target cells. CD8+ T-cells include most of the cytotoxic T-cells.

**References**

1. Knapp W. et. al. Leukocyte Typing IV, p342-343, (1989)
2. Parnes JR, Semin Immunol 6: 221-229 (1994)
3. Delon J. et al. Immunity 9(4): 467-73 (1998)
4. Akimoto H, et al. Immunology 95(2): 214-218 (1998)
5. Leahy DJ. Transplantation 45(4): 677-682 (1998)

**FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES**