



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Mouse anti-Uroplakin III, clone AU1 (Monoclonal)

Clone no. AU1

MONOSAN

---

Product name	Mouse anti-Uroplakin III, clone AU1 (Monoclonal)
Host	Mouse
Applications	IHC-P (1:50-1:100)
Species reactivity	human
Conjugate	-
Immunogen	AUM preparation from bovine urinary bladder
Isotype	IgG1
Clonality	Monoclonal
Clone number	AU1
Size	1 ml
Concentration	lot specific
Format	-
Storage buffer	Tissue culture supernatant with sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Uroplakin III, clone AU1 (Monoclonal)

Clone no. AU1

MONOSAN

**Additional info**

Uroplakins (UPs) are a family of transmembrane proteins (UPs Ia, Ib, II and III) that are specific differentiation products of urothelial cells. In non-neoplastic mammalian urothelium, UPs are expressed in the luminal surface plasmalemma of superficial (umbrella) cells, forming complexes of 16nm crystalline particles. UPIII is specific for tumors of urothelial origin and, when used in combination with other markers, can aid in the diagnosis of primary and metastatic tumors.

**References**

1. Moll et al. Am J Pathol 1995;147:1383-1397
2. Olsburgh et al. J Pathol 2003;199:41-49
3. Parket et al. 2003;27:1-10
4. Ohtsuka et al. 2006;97:1322-1326
5. -

**FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES**