



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Mouse anti-CA125 (Ovarian Cancer Antigen), clone OV185 (monoclonal)

Clone no. Ov185:1

MONXtra

Product name	Mouse anti-CA125 (Ovarian Cancer Antigen), clone OV185 (monoclonal)
Host	Mouse
Applications	IHC-P (1:100)
Species reactivity	human
Conjugate	-
Immunogen	A partially purified mucin fraction from a pool of cancer tissue from patients with epithelial ovarian cancer.
Isotype	IgG1
Clonality	Monoclonal
Clone number	Ov185:1
Size	1 ml
Concentration	Greater than or equal to 12.5 mg/L
Format	-
Storage buffer	Purified in PBS, 1% BSA and 15 mM sodium azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-CA125 (Ovarian Cancer Antigen), clone OV185 (monoclonal)

Clone no. Ov185:1

MONXtra

Additional info

CA125 antigen is usually associated with ovarian epithelial malignancies. Serum assays are widely used to detect this protein in the monitoring of ovarian cancers. CA125 antigen may also be detected by immunohistochemistry and expression has been found in neoplasms such as seminal vesicle carcinoma and anaplastic lymphoma. CA125 antigen is not found exclusively in malignant tumors. CA125 is also known as MUC16.

References

1. Weng D et al. Int. Journal of Cancer. 2011; 129:1990-2001
2. Gabriel M et al. Gynakol Geburtshilfliche Rundsch. 2000; 40(3-4):140-144
3. Gabriel M et al. Ginekol Pol. 1999; 70(11):819-823
4. Keen CE et al. Histopathology. 1999; 34:113-117
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES