



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Mouse anti-Utrophin, clone DRP3/20C5 (monoclonal)

Clone no. DRP3/20C5

MONXtra

---

Product name	Mouse anti-Utrophin, clone DRP3/20C5 (monoclonal)
Host	Mouse
Applications	IHC-fr
Species reactivity	human, rat, dog
Conjugate	-
Immunogen	Amino terminal domain
Isotype	IgG1
Clonality	Monoclonal
Clone number	DRP3/20C5
Size	2.5ml
Concentration	n/a
Format	-
Storage buffer	Lyophilized
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Mouse anti-Utrophin, clone DRP3/20C5 (monoclonal)

Clone no. DRP3/20C5

MONXtra

**Additional info**

The utrophin gene is located on chromosome 6. Analyte Specific Reagent. Analytical and performance characteristics are not established. Amino terminal domain of the human homolog of human dystrophin, utrophin (also known as dystrophin related protein or "DRP"). Also crossreacts with utrophin in sections of muscle from rat and dog. Other animals species have not been tested. The product 2 is a lyophilized tissue culture supernatant containing sodium azide as a preservative. The user is required to reconstitute the contents of the vial with the correct volume of sterile distilled water as indicated on the vial label.

**References**

1. -
2. -
3. -
4. -
5. -

**FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES**