



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Product datasheet PS210

MONOSAN[®]

Polyclonal Rabbit anti-Substance P 4-11

Clone no. -

MONOSAN

Product name	Polyclonal Rabbit anti-Substance P 4-11
Host	Rabbit
Applications	IHC-fr, IHC-P (1:400-1:800)
Species reactivity	Human, Rat
Conjugate	-
Immunogen	Synthetic Substance P 4-11, conjugated to thyroglobulin using CDI
Isotype	-
Clonality	Polyclonal
Clone number	-
Size	50 µl (lyoph.)
Concentration	n/a
Format	-
Storage buffer	Lyophilized; reconstitute in 100 µl dist. water
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Polyclonal Rabbit anti-Substance P 4-11

Clone no. -

MONOSAN

Additional info

Substance P occurs in nerve fibers of the central and peripheral nervous system and in endocrine cells of the gut. It stimulates smooth muscle contraction, gives rise to vasodilation and is involved in sensory functions. Substance P-containing tumors arising in the ileum are often associated with the carcinoid syndrome, characterized by flushing of the skin, diarrhea, broncho-constriction and sudden drops in blood pressure. Substance P is commonly found in the midgut carcinoids and some of the symptoms may be related to this peptide.
Absorption with 10-100 ug SP 1-4 per ml diluted antiserum abolishes the staining.
Positive control: frozen sections of rat colon.

References

1. -
2. -
3. -
4. -
5. -

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES