



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Product datasheet PS318

**MONOSAN**<sup>®</sup>

Polyclonal antibody against Granzyme B

Clone no. -

MONOSAN

---

Product name	Polyclonal antibody against Granzyme B
Host	Rabbit
Applications	IHC-P (1:10-1:50)
Species reactivity	human
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	-
Clonality	Polyclonal
Clone number	-
Size	1 ml
Concentration	n/a
Format	-
Storage buffer	Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Polyclonal antibody against Granzyme B

Clone no. -

MONOSAN

---

**Additional info**

Granzymes are serine proteases which are stored in specialized lytic granules of cytotoxic T lymphocytes and in natural killer cells. Anti-Granzyme B has been useful in diagnosing Natural killer/T cell lymphoma, as well as anaplastic large cell lymphoma. High percentages of cytotoxic T cells have been shown to be an unfavorable prognostic indicator in Hodgkin's disease.

**References**

1. Kummer JA, et al. Clin Exp Immunol. 1995; 100:164-72
2. Oudejans JJ, et al. Am J Pathol. 1996; 148:233-40
3. Oudejans JJ, et al. Blood. 1997; 89:1376-82
4. Kato N, et al. Am J Dermatopathol. 2003; 25:142-7
5. Liu J, et al. J Dermatol. 2003; 30:735-41

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES