



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Product datasheet PS333

**MONOSAN**<sup>®</sup>

Polyclonal antibody against Myeloperoxidase

Clone no. -

MONOSAN

---

Product name	Polyclonal antibody against Myeloperoxidase
Host	Rabbit
Applications	IHC-P (1:100-1:500)
Species reactivity	human
Conjugate	-
Immunogen	Unknown or proprietary to MONOSAN and/or its suppliers
Isotype	-
Clonality	Polyclonal
Clone number	-
Size	1 ml
Concentration	n/a
Format	-
Storage buffer	Tris Buffer, pH 7.3-7.7, containing 1% BSA and <0.1% Sodium Azide
Storage until expiry date	2-8°C

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES

Polyclonal antibody against Myeloperoxidase

Clone no. -

MONOSAN

---

**Additional info**

Anti-Myeloperoxidase detects granulocytes and monocytes in blood and precursors of granulocytes in the bone marrow. This antibody can detect myeloid cell populations of the bone marrow as well as in other sites.

**References**

1. Pinkus GS, et al. Mod Pathol. 1991; 4:733-41
2. Markoc F, et al. Tumori. 2010; 96:149-53
3. Alexiev BA, et al. Diagn Pathol. 2007; 31;2:42
4. Saravanan L, et al. Int J Lab Hematol. 2010; 32:132-6
5. Manaloor EJ, et al. Am J Clin Pathol. 2000; 113:814-22

FOR RESEARCH USE ONLY. NOT FOR USE IN DIAGNOSTIC PROCEDURES