



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Polyclonal Antibody to Human Kappa Light Chain

Description:

- Immunogen:** Kappa Light chain isolated from pooled urine of patients with Bence Jones proteinuria.
- Format:** This antibody has been pretitered and quality controlled to work on formalin-fixed paraffin-embedded as well as acetone fixed cryostat tissue sections. No further titration is required.
- Specificity:** This antibody reacts with free as well as bound kappa light chains. Contaminating antibodies have been removed by solid phase absorption.

Uses/Limitations:

Immunohistochemistry
For Research Use Only.
Do not use past expiration date.

Storage:

2-8° Centigrade.

Procedure:

We suggest an incubation period of 30-60 minutes at room temperature. However, depending upon the fixation conditions and the staining system employed, optimal incubation should be determined by the user.

Precautions:

Contains Sodium Azide as a preservative.
Do not pipette by mouth.