



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Bluing Reagent For Hematoxylin

Description: Bluing Reagent is designed to provide complete bluing of hematoxylin without damaging cell structure.

Results: Nuclei are blue.

Availability:

<u>Item #</u>	<u>Volume</u>
BRT999	1000ml
BRT500	500 ml

Procedure:

1. Stain with hematoxylin according to protocol.
2. Place in Bluing Reagent until nuclei are blue (1-2 minutes).
3. Rinse in DI water.
4. Dehydrate, clear, and mount in permanent mounting medium.