



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

Wash Solution Concentrate (50x)

Description: Wash Solution Concentrate (50x) is designed to provide optimal rinsing of routine or special stains between steps. In addition, this reagent can be used as a dilution buffer for most aqueous stains. This product is compatible with hand staining procedures and most automated systems.

Appearance: Transparent liquid, slightly yellow. Discard if reagent becomes cloudy.

Uses: Routine and Special Stain procedures.
Stain diluent.

Limitations: Do not use past expiration date.

Availability:

<u>REF</u>	<u>Volume</u>
WSC-1	4 x 200 ml
WSC999	1000 ml
WSC20000	20 liters

Storage: Store at 18-25°C. Product is stable for 24 months from the date of manufacture.

Procedure:

1. Add 200ml of Wash Solution Concentrate (50x) to 9.8 liters of DI/Distilled water and mix thoroughly at low speed to avoid “foaming” the reagent.
2. Mix fresh daily for optimal results.

Precautions: Do not pipette by mouth.
Avoid contact with skin and eyes.
Observe all federal, state and local environmental regulations regarding disposal.

Storage: 18° C



25° C



ScyTek Laboratories, Inc.
205 South 600 West
Logan, UT 84321
U.S.A.



EMERGO
EmergoEurope (31)(0) 70 345-8570
Molsnstraat 15
2513 BH Hague, The Netherlands