

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere Liefer- und Versandbedingungen

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

linkedin.com/company/szaboscandic in



Data Sheet (Cat.No.T1178)



Temozolomide

Chemical Properties

CAS No.: 85622-93-1

Formula: C6H6N6O2

Molecular Weight: 194.15

Appearance: no data available

Storage: store under nitrogen,keep away from direct sunlight

Powder: -20°C for 3 years | In solvent: -80°C for 1 year

Biological Description

| Description | Temozolomide (TMZ) is a DNA alkylating agent with blood-brain barrier permeability and oral activity. Temozolomide has antitumor activity and antiangiogenic activity, and also induces apoptosis and autophagy. Temozolomide is stable under acidic condition and hydrolyzes under neutral or slightly alkaline conditions. | | |
|---------------|--|--|--|
| Targets(IC50) | Apoptosis,DNA Alkylator/Crosslinker,DNA/RNA Synthesis,Autophagy | | |
| In vitro | METHODS: Melanoma cells SK-mel-28, MM200, IgR3, Mel-FH were treated with Temozolomide (0-500 μM) for 72 h. Cell viability was examined using MTT. RESULTS: The p53 status and MGMT expression levels were correlated with the sensitivity of Temozolomide. MM200 and IgR3 (expressing wild-type p53 and low MGMT levels) showed comparable sensitivity to Temozolomide, with IC50 values of 23 and 22 μM, respectively, whereas SK mel-28 and Mel-FH (mutant-type p53 and high MGMT level) were resistant with IC50 values >256 and >247 μΜ. [1] METHODS: Melanoma cells MM200 and IgR3 were treated with Temozolomide (100 μM) for 24-72 h. The cell cycle was examined by Flow Cytometry. RESULTS: Temozolomide induced G2/M cell cycle arrest in MM200 and IgR3 cells. [1] METHODS: Human glioma cells U118 were treated with Temozolomide (250-500 μM) for 3-48 h. The m5C level in DNA was measured. RESULTS: The response of U118 cells to Temozolomide depends on the concentration and time of the reaction. The amount of m5C in DNA increased significantly within a short period of time after Temozolomide treatment. m5C(R) reached the highest level after 24 h of treatment with 500 μM Temozolomide. [2] | | |
| In vivo | METHODS: To assay antitumor activity in vivo, Temozolomide (68 mg/kg by gavage) and AG-014699 (1 mg/kg by intraperitoneal injection) were administered intraperitoneally to CD1 nu/nu mice harboring medulloblastomas D425Med, D283Med, or D384Med once daily for five days. RESULTS: AG-014699 enhanced the efficacy of Temozolomide in an in vivo model of medulloblastoma. [3] METHODS: To assay antitumor activity in vivo, Temozolomide (0.9 mg/kg orally once daily) and Aldox (16 mg/kg intravenously once weekly) were administered to Foxn1 nude mice bearing human glioblastoma U87MG once daily for five weeks. RESULTS: Combined treatment with Temozolomide and AldoxAldo induced significant tumor volume suppression and increased survival. [4] | | |

Page 1 of 2 www.targetmol.com

| Cell Research | Cell lines exposed to TMZ (with or without 5-Aza or O6-BG pre-treatment) were grown in 24-well plates under standard culture conditions for 6 days. Cytotoxicity was determined using the sulphorhodamine-B (SRB) method. Briefly, the cells were fixed with 10% trichloroacetic acid for 20 min at 4°C then washed three times with water. After 24 hours, cells were stained for 30 min at room temperature with 0.4% SRB dissolved in 1% acetic acid and then washed three times with 1% acetic acid. The plates were air-dried and the dye solubilized with 300 ml/well of 10 mM Tris base (pH 10.5) for 10 min on a shaker. The optical density of each well was measured spectrophotometrically using a Titertek multiscan colorimeter at 492 nm [2]. | |
|-----------------|---|--|
| Animal Research | TZM was dissolved in dimethyl-sulfoxide (40 mg/mL), diluted in saline (5 mg/mL), and administered intraperitoneally on day 2 after tumor injection at 100 mg/kg or 200 mg/kg, doses commonly used for in vivo preclinical studies.15-17 Because cytotoxicity induced by TZM and PARP inhibitors can be improved by fractionated modality of treatment,9 in selected groups a total dose of 200 mg/kg TZM was divided in 2 doses of 100 mg/kg given on days 2 and 3. NU1025 was dissolved in polyethylene glycol-400 (40% in saline) and was injected intracranially at the maximal deliverable dose (1 mg/mouse, 0.03 mL) or, in selected groups, intraperitoneally (0.3 mL) on day 2 after tumor challenge, 1 hour before TZM administration. Control mice were injected with drug vehicles [4]. | |

Solubility Information

| Solubility | PBS: 5 mg/mL (25.75 mM),Sonication is recommended.(The compound is unstable in | | |
|------------|--|--|--|
| | PBS. Please use soon) | | |
| | DMSO: 28.83 mg/mL (200 mM) | | |
| | (< 1 mg/ml refers to the product slightly soluble or insoluble) | | |

Preparing Stock Solutions

| | 1mg | 5mg | 10mg |
|-------|-----------|------------|------------|
| 1 mM | 5.1507 mL | 25.7533 mL | 51.5066 mL |
| 5 mM | 1.0301 mL | 5.1507 mL | 10.3013 mL |
| 10 mM | 0.5151 mL | 2.5753 mL | 5.1507 mL |
| 50 mM | 0.103 mL | 0.5151 mL | 1.0301 mL |

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Wang Y, Wang X, Wang K, et al. Chronic stress accelerates glioblastoma progression via DRD2/ERK/ β -catenin axis and Dopamine/ERK/TH positive feedback loop. Journal of Experimental & Clinical Cancer Research. 2023, 42(1): 1-

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only. Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street, Wellesley Hills, MA 02481

Page 2 of 2 www.targetmol.com