



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com


www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic)

AGLME

Chemical Properties

CAS No. :	10236-44-9
Formula:	C ₁₁ H ₂₁ N ₃ O ₄
Molecular Weight:	259.3
Appearance:	no data available
Storage:	keep away from moisture Powder: -20°C for 3 years In solvent: -80°C for 1 year



Biological Description

Description	AGLME is used in a direct enzymatic assay for activated Hageman factor measuring the ability of Hageman factor to hydrolyze the cpd.
-------------	--

Preparing Stock Solutions

	1mg	5mg	10mg
1 mM	3.8565 mL	19.2827 mL	38.5654 mL
5 mM	0.7713 mL	3.8565 mL	7.7131 mL
10 mM	0.3857 mL	1.9283 mL	3.8565 mL
50 mM	0.0771 mL	0.3857 mL	0.7713 mL

Please select the appropriate solvent to prepare the stock solution, according to the solubility of the product in different solvents. Please use it as soon as possible.

Reference

Capet-Antonini FC, Tamenasse J. Kinetic studies on soluble and insoluble urokinases. Can J Biochem. 1975 Aug;53(8):890-4. PubMed PMID: 241473.

Inhibitor · Natural Compounds · Compound Libraries · Recombinant Proteins

This product is for Research Use Only · Not for Human or Veterinary or Therapeutic Use

Tel:781-999-4286 E_mail:info@targetmol.com Address:36 Washington Street,Wellesley Hills,MA 02481