



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com



www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 



Goat Anti-Fluorescein, Biotinylated BA-0601-.5

Product Images



Short Description

Fluorescein is not found in biological systems. This green fluorescent label can be directly visualized or used as a specific tag and detected with our high affinity, purified antibody reagents. Labeled probes can be used in applications such as *in situ* hybridization and blotting.

Additional Information

Unit Size	0.5 mg
Applications	In situ hybridization, Blotting Applications
Concentration	0.5 mg active conjugate/ml
Recommended Storage	2-8 °C
Solution	10 mM sodium phosphate, 0.15 M NaCl, pH 7.8, 0.08% sodium azide, 3 mg/ml bovine serum albumin.
Recommended Usage	The recommended concentration range for use is 2 - 10 µg/ml. If this biotinylated antibody is to be used in tissues, which may contain cross-reacting endogenous immunoglobulins, dilution of this biotinylated antibody may be made in buffers containing 2% normal serum from the same species as the tissue.
Detection Target	Fluorescein
Conjugate	Biotinylated
Host Species	Goat
Format	Concentrate

