



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 



Lectin Kit III, Biotinylated

BK-3000

Product Images



Additional Information

Unit Size	1 kit
Applications	Immunohistochemistry / Immunocytochemistry, Immunofluorescence, Blotting Applications, Elispot, ELISAs, Glycobiology
Recommended Usage	Reconstitute lyophilized biotinylated lectins by adding 0.5 ml water. The resulting solution will be: 10 mM HEPES, 0.15 M NaCl, 0.1 mM CaCl ₂ , 0.08% sodium azide, pH 7.5.
Conjugate	Biotinylated
Sugar Specificity	Galactose, N-Acetylgalactosamine, N-Acetylglucosamine

