



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 



Datura Stramonium Lectin (DSL), Unconjugated

L-1180-5

[Product Images](#)



Short Description

The carbohydrate binding site recognizes (β -1,4) linked N-acetylglucosamine oligomers, preferring chitobiose or chitotriose over a single N-acetylglucosamine residue. This lectin binds well in the acidic pH range but its affinity decreases above pH 8.0. DSL also binds well to N-acetyllactosamine and oligomers containing repeating N-acetyllactosamine sequences. A branched pentasaccharide including two N-acetyllactosamine disaccharides linked to mannose (β -1,6) and (β -1,2) was reported to be the most potent inhibitor of agglutination.

Additional Information

Unit Size	5 mg
Applications	Glycobiology, Mitogenic Stimulation
Recommended Usage	Though many buffers can be employed for reconstituting and diluting this lectin, 10 mM HEPES buffered saline, pH 8.5, 0.1 mM CaCl_2 is recommended. For preserving solutions stored at 4 °C, 0.08% sodium azide can be used. Aggregation may occur with time if stored at concentrations greater than 2 mg/ml.
Recommended Storage	2-8 °C; for long term storage, aliquots may be stored frozen or preserved with 0.08% sodium azide in the recommended buffer and stored at 2-8 °C
Conjugate	Unconjugated
Sugar Specificity	[GlcNAc]1-3, N-Acetylglucosamine

