



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 



ZYMO RESEARCH

The Beauty of Science is to Make Things Simple

INSTRUCTION MANUAL

ZR 1 kb DNA Marker™

Catalog Nos. M5003 & M5006

Contents

I. General Information.....	1
Description	
Storage Conditions	

For Research Use Only

ZYMO RESEARCH CORP.

Phone: (949) 679-1190 ▪ Toll Free: (888) 882-9682 ▪ Fax: (949) 266-9452 ▪ info@zymoresearch.com ▪ www.zymoresearch.com

I. General Information

Description

Zymo Research's **DNA Molecular Weight Markers™** are suitable for use as molecular weight standards for agarose and polyacrylamide gel electrophoresis. The 1 kb DNA marker consists of plasmid DNA that has been digested with endonucleases to yield 10 defined fragments that range in size from 0.5 - 10 kb. The intensity of the 3 kb band is increased as an identifier. All fragments have 4-base, 5' overhangs.

Effective Size Range: 500 bp to 10 kb

Source:

The double-stranded DNA is digested to completion with appropriate restriction enzymes, phenol extracted and equilibrated to 10 mM Tris-HCl (pH 8.0) and 1 mM EDTA.

Concentration:

500 µg/ml

Recommended Load: 0.5 µg

Storage Conditions

Storage Conditions:

10 mM Tris-HCl

1 mM EDTA

Storage Temperature:

4°C

For long term storage (>30 days), store at -20°C.