

# Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



## Lieferung & Zahlungsart

siehe unsere Liefer- und Versandbedingungen

# Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

#### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

linkedin.com/company/szaboscandic in





# Anti-Outer membrane protein [CB-NM.2] Bulk Size Ab03266-23.0-BT

This chimeric rabbit antibody was made using the variable domain sequences of the original Mouse IgG2b format for improved compatibility with existing reagents assays and techniques.

Isotype and Format: Rabbit IgG, Kappa

Clone Number: CB-NM.2

Alternative Name(s) of Target: class 3 outer membrane protein; 15F9; NM-2

**UniProt Accession Number of Target Protein:** 

Published Application(s): WB, ELISA

Published Species Reactivity: Neisseria meningitidis

**Immunogen:** The original antibody was generated by immunizing mice with serotype antigen preparations from case isolate B385 strain of Neisseria meningitidis, from Havana.

**Specificity:** This antibody is specific for a five amino acid region (S/T)VETG located in the main variable region (VR) VR1 which was present in all N. meningitidis strains. It recognizes class 3 outer membrane protein of Neisseria meningitidis. This antibody can recognize fifty six strains belonging to seven serogroups: A, B, 29E, L, X, Y, Z, and five serotypes: 1, 4, 5, 12, 13 but did not recognize the serotypes 2, 6, 9, 11, 14, 15, 18 and 19.

**Application Notes:** The binding reactivity of this antibody against fifty six strains of Neisseria meningitidis belonging to seven serogroups: A, B, 29E, L, X, Y, Z, and five serotypes: 1, 4, 5, 12, 13 was assessed using whole cell ELISA. This antibody can also be used in the detection of these serogroups and serotypes using western blot (Guillen et al., 1997).

**Antibody First Published in:** Guillen et al. Molecular analysis Neisseria meningitidis class 3 outer membrane protein in strains recognized by the monoclonal antibody CB-Nm.2. Biotechnologia Applicada 1997; Vol -14:23-30. PMID:

**Note on publication:** This paper analyzes the reactivity of this antibody with various serogroups and serotypes of N. meningitidis and determines the binding epitope of this antibody.

### **Product Form**

**Size:** 1 mg Purified antibody in bulk size.

**Purification:** 

Protein A affinity purified **Supplied In:** PBS only.

**Storage Recommendation:** Store at 4°C for up to 3 months. Note, this antibody is provided without added preservatives, it is therefore recommed this antibody be handled under sterile conditions. For longer storage, aliquot and store at -20°C.

**Concentration:** 1 mg/ml.

Important note – This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals.