



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

## Anti-P32 [AbAb-LSD-1] Bulk Size Ab03976-23.0-BT

**Isotype and Format:** Rabbit IgG, Kappa

**Clone Number:** AbAb-LSD-1

**Alternative Name(s) of Target:** LSD; Lumpy skin disease

**UniProt Accession Number of Target Protein:** A0A6M8DVW4

**Published Application(s):** ELISA

**Published Species Reactivity:** bovine nodular skin disease virus

**Immunogen:** The original antibody was generated by immunizing female Balb/c mice with LSDV-P32 protein antigen.

**Specificity:** The antibody is specific for P32.

**Application Notes:** The specificity of the antibody for P32 was confirmed by ELISA analysis (CN112111006B).

**Antibody First Published in:** [PMID:](#)

**Note on publication:**

### Product Form

**Size:** 1 mg Purified antibody in bulk size.

**Purification:** Protein A affinity purified

**Supplied In:** PBS only.

**Storage Recommendation:** Store at 4°C for up to 3 months. Note, this antibody is provided without added preservatives, it is therefore recommended this antibody be handled under sterile conditions. For longer storage, aliquot and store at -20°C.

**Concentration:** 1 mg/ml.

Important note - This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals.