



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

## Anti-VWFpp [JA9] Standard Size Ab04116-30.11

This is an scFv fragment with a His tag.

This is a reformatted human scFv antibody, based on the original human scFv format, created for improved compatibility with existing reagents, assays and techniques.

**Isotype and Format:** scFv fragment (His), ScFv

**Clone Number:** JA9

**Alternative Name(s) of Target:** von Willebrand factor propeptide

**UniProt Accession Number of Target Protein:**

**Published Application(s):** ELISA

**Published Species Reactivity:** Human

**Immunogen:** The original antibody was generated by selecting antibody fragments using phage display technology, where two phage display libraries were amplified and used to screen for binding to von Willebrand factor propeptide (VWFpp) displayed on yeast cells.

**Specificity:** This antibody is specific for VWF propeptide (VWFpp) and not to the yeast used for the expression of the propeptides.

**Application Notes:** This antibody's original format (scFv) was used in a sandwich ELISA to determine the VWF propeptide levels in plasma and showed a higher binding affinity for VWFpp in plasma at especially lower plasma concentrations (Meiring et al., 2018; PMID: 29992156).

**Antibody First Published in:** Meiring et al. The Use of Phage Display and Yeast Based Expression System for the Development of a Von Willebrand Factor Propeptide Assay: Development of a Von Willebrand Factor Propeptide Assay Biomed Res Int. 2018 May 24;2018:6232091. doi: 10.1155/2018/6232091 [PMID:29992156](#)

**Note on publication:** The original publication describes the development of a cost-effective assay for the diagnosis of von Willebrand disease, focusing on the detection of von Willebrand factor propeptide using phage display and yeast-based expression systems.

## Product Form

**Size:** 100 µg Purified antibody.

**Purification:** Purified by Immobilized Metal Affinity Chromatography

**Supplied In:** PBS with 0.02% Proclin 300.

**Storage Recommendation:** Store at 4°C for up to 3 months. For longer storage, aliquot and store at -

20°C.

**Concentration:** 1 mg/ml.

Important note - This product is for research use only. It is not intended for use in therapeutic or diagnostic procedures for humans or animals.