

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere Liefer- und Versandbedingungen

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

linkedin.com/company/szaboscandic in





IGH@ FISH Probe

Catalog #: FA0341 規格:[200 uL]

List All

Specification	
Product Description:	Made to order FISH probes for identification of gene amplification using Fluorescent In Situ Hybridization Technique. (<u>Technology</u>)
Supplied Product:	DAPI Counterstain (1500 ng/mL) 250 uL
Storage Instruction:	Store at 4°C in the dark.
Origin:	Human
Source:	Genomic DNA
Notice:	We strongly recommend the customer to use FFPE FISH PreTreatment Kit 1 (Catalog #: KA2375 or KA2691) for the pretreatment of Formalin-Fixed Paraffin-Embedded (FFPE) tissue sections.
Regulation Status:	For research use only (RUO)
A!:4:	
Applications	01. II. II. II. II. II. II. II. II. II. I
Fluorescent in	Situ Hybridization (Cell)
Gene Information	on
Entrez GenelD:	3492
Gene Name:	IGH
Gene Name: Gene Alias:	IGH IGH,IGH.1@,IGHDY1,MGC72071,MGC88774
Gene Alias:	
Gene Alias:	IGH,IGH.1@,IGHDY1,MGC72071,MGC88774 immunoglobulin heavy locus

which is alternatively spliced to encode either a mu or a delta heavy chain. Mature B cells in the lymph nodes undergo switch recombination, so that the V-D-J gene is brought in proximity to one of the IGHG, IGHA, or IGHE genes and each cell expresses either the gamma, alpha, or

Application Image

Fluorescent In Situ Hybridization (Cell)

2017/2/16

epsilon heavy chain. Recombination of many different V segments with several J segments provides a wide range of antigen recognition. Additional diversity is attained by junctional diversity, resulting from the random additional of nucleotides by terminal deoxynucleotidyltransferase, and by somatic hypermutation, which occurs during B cell maturation in the spleen and lymph nodes. Several V, D, J, and C segments are known to be incapable of encoding a protein and are considered pseudogenes. [provided by RefSeq

Other -Designations:

Related Disease

Chromosome Aberrations

服務條款 | 隱私權政策 | 著作及商標 | 網站地圖

©2017 亞諾法生技股份有限公司 Abnova Corporation. 版權所有.

Page 2 of 2 2017/2/16