



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 

identify

-

preserve

-

process

AlteRead™

2D barcode reader set-up packs

Store & track samples easily with all the kit supplied in one cost saving package

AlteRead™ Large & SBS rack set-up pack



Compatible with

81, 100, 192 Large racks & all SBS formats;
12, 24, 48, 96 & 384 racks

AlteRead™ SBS rack set-up pack



Compatible with

24, 48 & 96 SBS rack formats

Each AlteRead set-up pack includes:

- AlteRead™ rack 2D barcode camera reader
- AlteRead™ Single tube 2D barcode camera reader
- Windows laptop with AlteRead software pre-loaded
- 2D barcoded Large format or SBS AlteTubes
- 2 cases of tubes, racks & caps in any packaging format
- AlteTube™ Large format or SBS tube selector

AlteRead™ Single 2D barcode camera reader



Versatile

Easy to use software included & display results in any application


Small and portable

Powered by a single USB cable

Fast decoding

Decode a single tube 2D barcode in less than 1 second

AlteRead™ rack 2D barcode camera readers



Instant results

Capture, read, decode and export data for a full rack in less than 2 seconds

Easy to use

Includes pre-calibrated intuitive AlteRead software, ideal for multi-user labs

Small footprint

Uses little space on the bench or when integrated on the deck of robotic systems

Cryo Cover™

Reduces condensation for uninterrupted reading of cold racks

LIMS integration

Seamlessly transfer AlteRead™ export data into databases

