

# Produktinformation



Forschungsprodukte & Biochemikalien
Zellkultur & Verbrauchsmaterial
Diagnostik & molekulare Diagnostik
Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



Lieferung & Zahlungsart siehe unsere Liefer- und Versandbedingungen

# Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

# SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien T. +43(0)1 489 3961-0 F. +43(0)1 489 3961-7 <u>mail@szabo-scandic.com</u> www.szabo-scandic.com





# Interleukin-17A/F Heterodimer, rat recombinant (rrIL-17A/F)

Catalog No:	97291
Lot No:	XXXXX
Source:	E. coli
Synonyms:	IL17A/F, IL17 A/F, IL-17A/F, IL-17 A/F, IL17AF, IL-17 AF, Interleukin-17 A/F, Interleukin-17 AF

# Background

Rat IL-17A/F is a glycoprotein which is secreted as a disulfide-linked heterodimer. IL-17A/F consists of two proteins of the IL-17 family, IL-17A and IL17F. Proteins of the 6 homodimeric IL17 family show a cysteine knot motif that contains two disulfide-bonds. Rat IL17A is produced as a 155 a.a precursor that includes a 23 amino acids signal sequence and a 132 amino acid chain that includes an N-linked glycosylation site. Rat IL17F is produced as a 153 amino acid precursor with a 20 amino acid signal sequence and a 133 amino acid region. Similar to IL17A, IL17F also has an N-linked glycosylation site. Both proteins (IL17A & IL17F) share 50% amino acid sequence identity. Human IL17A & IL17F show approximately 60% homology in their amino acid sequence to mouse IL-17A and IL-17F. Interleukin-17A/F and IL17A, IL17F homodimers are manufactured by activted CD4+ T cells, called Th17. IL-23 causes Th17 lymphocytes to manufacture IL-17A/F. IL17RA and IL17RC form a heterodimer for the binding of IL17A and IL17F. IL-17A/F binds IL-17RA. Interleukin-17A/F induces chemokine production and airway neutrophilia with intermediate potency between IL17A (most potent) and IL17F (least potent).

## Description

Interleukin-17A/F rat recombinant produced in *E. coli* is a heterodimeric, non-glycosylated polypeptide chain containing 1 monomeric subunit of each IL-17A and IL-17F. The dimer contains 269 amino acids and having a total molecular mass of 30.7 kDa. IL-17A/F is purified by proprietary chromatographic techniques.

## **Physical Appearance**

Sterile filtered white lyophilized (freeze-dried) powder.

# Formulation

Lyophilized from a concentrated (1 mg/ml) solution containing no additives.

## Solubility

It is recommended to reconstitute the lyophilized Interleukin-17A/F in sterile water not more than 1 mg/ml, which can then be further diluted to other aqueous solutions.

## Stability

Lyophilized IL-17A/F, although stable at room temperature for 3 weeks, should be stored desiccated below -18°C. Upon reconstitution Rat IL-17A/F should be stored at 4°C between 2-7 days and for future use below -18°C. For long term storage it is recommended to add a carrier protein (0.1% HSA or BSA). Please prevent freeze-thaw cycles.

## Purity

Greater than 98.0% as determined by (a) HPLC and (b) SDS-PAGE.

## **Amino Acid Sequence**

IL-17A:	MAVLIPQSSV	CPNAEANNFL	QNVKVNLKVL	NSLSSKASSR	RPSDYLNRST	SPWTLSRNED	PDRYPSVIWE
	AQCRHQRCVN	AEGKLDHHMN	SVLIQQEILV	LKREPEKCPF	TFRVEKMLVG	VGCTCVSSIV	RHAS
IL-17B:	MARRNPKVGL	SALQKAGNCP	PLEDNSVRVD	IRIFNQNQGI	SVPRDFQNRS	SSPWDYNITR	DPDRFPSEIA
	EAQCRHSGCI	NAQGQEDGSM	NSVPIQQEIL	VLRREPQGCS	NSFRLEKMLI	KVGCTCVTPI	VHHAA

CONTACT US TODAY • BIOMOL GmbH • Kieler Str. 303A • 22525 Hamburg • Germany • info@biomol.com • www.biomol.com Fon: +49 (0)40-853 260 0 • Fax: +49 (0)40-853 260 22 • TOLL FREE IN GERMANY: Fon: 0800-246 66 51 • Fax: 0800-246 66 52





# **Biological Activity**

The ED50 as determined by its ability to induce IL-6 production from NIH-3T3 cells is 15.3 – 23 ng/ml.

# Usage

This product is offered by Biomol for research purposes only. Not for diagnostic purposes or human use. It may not be resold or used to manufacture commercial products without written approval of Biomol GmbH.