



# SZABO SCANDIC

Part of Europa Biosite

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!  
See the following pages for more information!



### Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

### Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

[mail@szabo-scandic.com](mailto:mail@szabo-scandic.com)

[www.szabo-scandic.com](http://www.szabo-scandic.com)

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 



6042 Cornerstone Court W, Ste B  
San Diego, CA 92121  
**Tel:** 1.858.829.3082  
**Fax:** 1.858.481.8694  
**Email:** [info@bpsbioscience.com](mailto:info@bpsbioscience.com)

## Thaw Medium 2 Catalog #: 60184-2

### DESCRIPTION:

This liquid Cell Thawing Medium is optimized for thawing and plating select BPS Bioscience cell lines, including

- PD-1/NFAT Reporter-Jurkat cell #60535
- NFAT Reporter (Luc) – Jurkat Cell Line #60621
- NF-κB-Luciferase Reporter (Luc) - Jurkat Cell Line #60651
- ADCC Bioassay Effector Cell F variant (Low Affinity) #60540
- ADCC Bioassay Effector Cell V variant (High Affinity) #60541
- GITR / NF-κB-Luciferase Reporter (Luc) - Jurkat Cell Line #60546

*Note: once cells are well-established, switch to the appropriate cell growth medium.*

### SIZE:

500 ml. Sterile filtered (0.2 μm) and Mycoplasma-free

### QUALITY:

Each lot is tested for growth and functionality of luciferase reporter-Jurkat cells. All lots are validated for low endotoxin levels (<1 EU/ml).

### STORAGE AND STABILITY:

This medium is shipped on gel packs. Upon receipt, the medium should be stored in the dark at 4°C. If frozen during shipping, please move to 4°C for storage. Protect from light. If the medium is warmed prior to use, do not exceed 37°C. When stored in the dark at 4°C, the product is stable until 6 months from date of receipt.

### THAWING METHOD:

It is recommended to quickly thaw the frozen cells from liquid nitrogen in a 37°C water-bath, then transfer the entire contents of the vial to a tube containing 10 ml of Cell Thawing medium. Spin down the cells, remove supernatant and resuspend the cells in pre-warmed Cell Thawing Medium. Transfer the resuspended cells to a T25 flask and incubate at 37°C in a 5% CO<sub>2</sub> incubator. Cells should be split before they reach 3x10<sup>6</sup> cells/ml.

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.829.3082** Fax **1.858.481.8694**

Or you can Email us at: [info@bpsbioscience.com](mailto:info@bpsbioscience.com)

Please visit our website at: [www.bpsbioscience.com](http://www.bpsbioscience.com)



6042 Cornerstone Court W, Ste B  
San Diego, CA 92121  
**Tel:** 1.858.829.3082  
**Fax:** 1.858.481.8694  
**Email:** [info@bpsbioscience.com](mailto:info@bpsbioscience.com)

## RELATED PRODUCTS:

<b>Product</b>	<b>Cat. #</b>	<b>Size</b>
PD-1 Cell Growth Medium	#60189	500 ml
PD-L1 Cell Growth Medium	#60188	500 ml
PD-1/NFAT Reporter-Jurkat cell	#60535	2 vials
NFAT Reporter (Luc) – Jurkat Cell Line	#60621	2 vials
NF-κB Luciferase Reporter (Luc) - Jurkat Cell Line	#60651	2 vials
ADCC Bioassay Effector Cell F variant (Low Affinity)	#60540	2 vials
ADCC Bioassay Effector Cell V variant (High Affinity)	#60541	2 vials
Thaw Medium 1	#60187	100 ml
Thaw Medium 3	#60186	100 ml
Thaw Medium 4	#60181	100 ml
Thaw Medium 5	#60182	100 ml
Thaw Medium 6	#60183	100 ml

OUR PRODUCTS ARE FOR RESEARCH USE ONLY. NOT FOR DIAGNOSTIC OR THERAPEUTIC USE.

To place your order, please contact us by Phone **1.858.829.3082** Fax **1.858.481.8694**

Or you can Email us at: [\*\*info@bpsbioscience.com\*\*](mailto:info@bpsbioscience.com)

Please visit our website at: [\*\*www.bpsbioscience.com\*\*](http://www.bpsbioscience.com)