

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere Liefer- und Versandbedingungen

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

linkedin.com/company/szaboscandic in



PRODUCT INFORMATION



Oleic Acid-biotin

Item No. 10010624

Formal Name: 9Z-octadecnoyl-N'-biotinoyl-1,5-

diaminopentane


MF: $C_{28}H_{50}N_4O_3S$

FW: 522.8 **Purity:** ≥95%

Supplied as: A solution in ethanol

-20°C Storage: Stability: ≥1 year

Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.



Laboratory Procedures

Oleic acid-biotin is supplied as a solution in ethanol. To change the solvent, simply evaporate the ethanol under a gentle stream of nitrogen and immediately add the solvent of choice. Solvents such as DMSO and dimethyl formamide purged with an inert gas can be used. The solubility of oleic acid-biotin in these solvents is approximately 50 mg/ml.

Further dilutions of the stock solution into aqueous buffers or isotonic saline should be made prior to performing biological experiments. Ensure that the residual amount of organic solvent is insignificant, since organic solvents may have physiological effects at low concentrations. If an organic solvent-free solution of oleic acid-biotin is needed, it can be prepared by evaporating the ethanol and directly dissolving the neat oil in aqueous buffers. The solubility of oleic acid-biotin in PBS (pH 7.2) is approximately 0.05 mg/ml. For greater aqueous solubility, linoleic acid-biotin can be directly dissolved in 0.1 M Tris-HCl (pH 8.5) (solubility of approximately 0.15 mg/ml) and then diluted with PBS (pH 7.2) to achieve the desired concentration or pH. We do not recommend storing the aqueous solution for more than one day. We do not recommend storing the aqueous solution for more than one day.

Description

Oleic acid is a monounsaturated fatty acid and is one of the major components of membrane phospholipids. Oleic acid contributes about 17% of the total fatty acids esterified to phosphatidylcholine, the major phospholipid class in porcine platelets. [3964] Oleic acid-biotin was designed to allow oleic acid to be detected in complexes with protein binding partners such as fatty acid binding proteins (FABPs). This it has the potential to serve as a tool to be used in the general elucidation of the signaling and transport of free oleic acid.

Reference

1. Wahle, K.W. and Peacock, L.I. Effects of isomeric cis and trans eighteen carbon monounsaturated fatty acids on porcine platelet function. Biochim. Biophys. Acta 1301(1-2), 141-149 (1996).

WARNING
THIS PRODUCT IS FOR RESEARCH ONLY - NOT FOR HUMAN OR VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.

This material should be considered hazardous until further information becomes available. Do not ingest, inhale, get in eyes, on skin, or on clothing. Wash thoroughly after handling. Before use, the user must review the complete Safety Data Sheet, which has been sent via email to your institution.

WARRANTY AND LIMITATION OF REMEDY

subject to Cayman's Terms and Conditions. Complete Terms and Conditions including Warranty and Limitation of Liability information can be found on our website

Copyright Cayman Chemical Company, 12/10/2021

CAYMAN CHEMICAL

1180 EAST ELLSWORTH RD ANN ARBOR, MI 48108 · USA PHONE: [800] 364-9897

[734] 971-3335

FAX: [734] 971-3640 CUSTSERV@CAYMANCHEM.COM WWW.**CAYMANCHEM**.COM