

# Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



# Lieferung & Zahlungsart

siehe unsere Liefer- und Versandbedingungen

# Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

# SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

linkedin.com/company/szaboscandic in



# PRODUCT INFORMATION



# **Cho-Arg (trifluoroacetate salt)**

Item No. 25943

CAS Registry No.: 1609010-56-1

Formal Name: cholest-5-en-3-ol (3β)-3-[6-[[(2S)-2-

amino-5-[(aminoiminomethyl)amino]-1-oxopentyl]amino]hexanoate],

2,2,2-trifluoroacetate (1:2) C<sub>39</sub>H<sub>71</sub>N<sub>5</sub>O<sub>3</sub> • 2CF<sub>3</sub>COO

MF: FW: 884.1


Supplied as: A solution in ethanol

≥95%

Storage: -20°C Stability: ≥1 year

**Purity:** 

Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.



#### **Laboratory Procedures**

Cho-Arg (trifluoroacetate salt) is supplied as a solution in ethanol. To change the solvent, simply evaporate the ethanol under a gentle stream of nitrogen and immediately add the solvent of choice. Solvents such as DMSO and dimethyl formamide purged with an inert gas can be used. The solubility of Cho-Arg (trifluoroacetate salt) in these solvents is approximately 50 mg/ml.

Cho-Arg (trifluoroacetate salt) is sparingly soluble in aqueous buffers. For maximum solubility in aqueous buffers, the ethanolic solution of Cho-Arg (trifluoroacetate salt) should be diluted with the aqueous buffer of choice. Cho-Arg (trifluoroacetate salt) has a solubility of approximately 0.25 mg/ml in a 1:3 solution of ethanol:PBS (pH 7.2) using this method.

### Description

Cho-Arg is a steroid-based cationic lipid that contains a cholesterol skeleton coupled to an L-arginine head group. It forms a complex with plasmid DNA and decreases plasmid DNA migration in an agarose-gel retardant assay at charge ratios greater than or equal to 4. Cho-Arg facilitates transfection of plasmid DNA into H1299 and HeLa cells in the presence or absence of fetal bovine serum, an effect that is reversed by the lipid raft-mediated endocytosis inhibitor methyl-β-cyclodextrin (Item No. 21633) and the caveolae-mediated endocytosis inhibitor genistein (Item No. 10005167). It is cytotoxic to H1299 cells (IC<sub>50</sub> = 88.5 μg/ml).

#### Reference

1. Sheng, R., Wang, Z., Luo, T., et al. Skeleton-controlled pDNA delivery of renewable steroid-based cationic lipids, the endocytosis pathway analysis and intracellular localization. Int. J. Mol. Sci. 19(2), e369 (2018).

WARNING
THIS PRODUCT IS FOR RESEARCH ONLY - NOT FOR HUMAN OR VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.

This material should be considered hazardous until further information becomes available. Do not ingest, inhale, get in eyes, on skin, or on clothing. Wash thoroughly after handling. Before use, the user must review the complete Safety Data Sheet, which has been sent via email to your institution.

## WARRANTY AND LIMITATION OF REMEDY

subject to Cayman's Terms and Conditions. Complete Terms and Conditions including Warranty and Limitation of Liability information can be found on our website

Copyright Cayman Chemical Company, 10/26/2018

### **CAYMAN CHEMICAL**

1180 EAST ELLSWORTH RD ANN ARBOR, MI 48108 · USA PHONE: [800] 364-9897

[734] 971-3335

FAX: [734] 971-3640 CUSTSERV@CAYMANCHEM.COM WWW.**CAYMANCHEM**.COM