

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere Liefer- und Versandbedingungen

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

linkedin.com/company/szaboscandic in



PRODUCT INFORMATION



Hexacosanoic Acid methyl ester-d₂

Item No. 28085

Formal Name: methyl hexacosanoate-26,26,26-d₃

C26:0 methyl ester-d₃, Ceric Acid methyl ester-d₃, Synonyms:

Cerotic Acid methyl ester-d₃

MF: $C_{27}H_{51}D_3O_2$ FW: 413.7


Chemical Purity: ≥98% (Hexacosanoic Acid methyl ester)

Deuterium

Incorporation: \geq 99% deuterated forms (d₁-d₃); \leq 1% d₀

Supplied as: A solid Storage: -20°C Stability: ≥2 years

Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.



Laboratory Procedures

Hexacosanoic acid methyl ester-d3 is intended for use as an internal standard for the quantification of hexacosanoic acid methyl ester (Item No. 26725) by GC- or LC-MS. The accuracy of the sample weight in this vial is between 5% over and 2% under the amount shown on the vial. If better precision is required, the deuterated standard should be quantitated against a more precisely weighed unlabeled standard by constructing a standard curve of peak intensity ratios (deuterated versus unlabeled).

Hexacosanoic acid methyl ester-d₃ is supplied as a solid. A stock solution may be made by dissolving the hexacosanoic acid methyl ester-d₃ in the solvent of choice, which should be purged with an inert gas. Hexacosanoic acid methyl ester-d3 is soluble in the organic solvent ethanol at a concentration of approximately 0.5 mg/ml.

Description

Hexacosanoic acid methyl ester is a fatty acid methyl ester. It has been used as a standard for the detection of hexacosanoic acid (Item No. 13354) in beeswax.¹

Reference

1. Du, J., Li, C., and Qiu, F.-P. GC-MS analysis of fatty acids in three kinds of natural beeswax. Shipin Keji 39(8), 268-271 (2014).

WARNING
THIS PRODUCT IS FOR RESEARCH ONLY - NOT FOR HUMAN OR VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.

This material should be considered hazardous until further information becomes available. Do not ingest, inhale, get in eyes, on skin, or on clothing. Wash thoroughly after handling. Before use, the user must review the complete Safety Data Sheet, which has been sent via email to your institution.

WARRANTY AND LIMITATION OF REMEDY

subject to Cayman's Terms and Conditions. Complete Terms and Conditions including Warranty and Limitation of Liability information can be found on our website

Copyright Cayman Chemical Company, 02/06/2020

CAYMAN CHEMICAL

1180 EAST ELLSWORTH RD ANN ARBOR, MI 48108 · USA PHONE: [800] 364-9897

[734] 971-3335

FAX: [734] 971-3640 CUSTSERV@CAYMANCHEM.COM WWW.**CAYMANCHEM**.COM