



SZABO SCANDIC

Part of Europa Biosite

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten!
See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere [Liefer- und Versandbedingungen](#)

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

[linkedin.com/company/szaboscandic](https://www.linkedin.com/company/szaboscandic) 


PRODUCT INFORMATION



Guanosine 5'-[γ-thio]triphosphate (lithium salt)

Item No. 35098

CAS Registry No.: 94825-44-2
Formal Name: guanosine 5'-(trihydrogen diphosphate), P'-anhydride with phosphorothioic acid, tetralithium salt
Synonyms: GTPγS, Guanosine 5'-O-(3-thio)triphosphate, Guanosine 5'-trihydrogen diphosphate monoanhydride, Phosphorothioic Acid
MF: C₁₀H₁₆N₅O₁₃P₃S • 4Li
FW: 567.0
Purity: ≥75%
UV/Vis.: λ_{max}: 253 nm
Supplied as: A solid
Storage: -20°C
Stability: ≥4 years



Special Conditions: Unstable in solution. Use immediately

Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

Laboratory Procedures

Guanosine 5'-[γ-thio]triphosphate (GTPγS) (lithium salt) is supplied as a solid. A stock solution may be made by dissolving the GTPγS (lithium salt) in water. The solubility of GTPγS (lithium salt) in water is approximately 75 mg/ml. We do not recommend storing the aqueous solution.

Description

GTPγS is a hydrolysis-resistant analog of GTP (Item No. 16060) that binds to and activates G proteins in the presence of magnesium.¹ GTPγS has been used to obtain crystal structures of GTPases in their active states, as well as measure the stability of GTPases in the presence of ligands.^{2,3} A radiolabeled form of GTPγS, [³⁵S]GTPγS, is commonly used in binding assays to determine agonism or antagonism at G protein-coupled receptors (GPCRs).^{1,4}

References

- Harrison, C. and Traynor, J.R. The [³⁵S]GTPγS binding assay: Approaches and applications in pharmacology. *Life Sci.* **74(4)**, 489-508 (2003).
- Ihara, K., Muraguchi, S., Kato, M., et al. Crystal structure of human RhoA in a dominantly active form complexed with a GTP analogue. *J. Biol. Chem.* **273(16)**, 9656-9666 (1998).
- Lito, P., Solomon, M., Li, L.-S., et al. Allele-specific inhibitors inactivate mutant KRAS G12C by a trapping mechanism. *Science* **351(6273)**, 604-608 (2016).
- Strange, P.G. Use of the GTPγS ([³⁵S]GTPγS and Eu-GTPγS) binding assay for analysis of ligand potency and efficacy at G protein-coupled receptors. *Br. J. Pharmacol.* **161(6)**, 1238-1249 (2010).

WARNING

THIS PRODUCT IS FOR RESEARCH ONLY - NOT FOR HUMAN OR VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.

SAFETY DATA

This material should be considered hazardous until further information becomes available. Do not ingest, inhale, get in eyes, on skin, or on clothing. Wash thoroughly after handling. Before use, the user must review the complete Safety Data Sheet, which has been sent via email to your institution.

WARRANTY AND LIMITATION OF REMEDY

Buyer agrees to purchase the material subject to Cayman's Terms and Conditions. Complete Terms and Conditions including Warranty and Limitation of Liability information can be found on our website.

Copyright Cayman Chemical Company, 12/01/2022

CAYMAN CHEMICAL

1180 EAST ELLSWORTH RD
ANN ARBOR, MI 48108 · USA

PHONE: [800] 364-9897

[734] 971-3335

FAX: [734] 971-3640

CUSTSERV@CAYMANCHEM.COM
WWW.CAYMANCHEM.COM