

## Produktinformation



Forschungsprodukte & Biochemikalien
Zellkultur & Verbrauchsmaterial
Diagnostik & molekulare Diagnostik
Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



Lieferung & Zahlungsart siehe unsere Liefer- und Versandbedingungen

## Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien T. +43(0)1 489 3961-0 F. +43(0)1 489 3961-7 <u>mail@szabo-scandic.com</u> www.szabo-scandic.com

# **PRODUCT** INFORMATION



## **Stearidonoyl Glycine**

Item No. 9000327

Formal Name:	2-(6Z,9Z,12Z,15Z)-octadeca- tetraenamidoacetic acid	
MF:	C <sub>20</sub> H <sub>31</sub> NO <sub>3</sub>	H
FW:	333.5	
Purity:	≥98%	
Supplied as:	A solution in ethanol	$\langle - \frown - \frown \rangle^0$
Storage:	-20°C	
Stability:	≥2 years	

Information represents the product specifications. Batch specific analytical results are provided on each certificate of analysis.

#### Laboratory Procedures

Stearidonoyl glycine is supplied as a solution in ethanol. To change the solvent, simply evaporate the ethanol under a gentle stream of nitrogen and immediately add the solvent of choice. Solvents such as DMSO and dimethyl formamide purged with an inert gas can be used. The solubility of stearidonoyl glycine in these solvents is approximately 20 and 30 mg/ml, respectively.

Stearidonoyl glycine is sparingly soluble in aqueous buffers. For maximum solubility in aqueous buffers, the ethanolic solution of stearidonoyl glycine should be diluted with the aqueous buffer of choice. The solubility of stearidonoyl glycine in PBS (pH 7.2) is approximately 1 mg/ml. We do not recommend storing the aqueous solution for more than one day.

#### Description

The  $\omega$ -3 polyunsaturated fatty acids (PUFAs) found in fish oils provide cardiovascular benefits. Stearidonic acid is an  $\omega$ -3 C18:4 PUFA. It is an intermediate in the metabolic pathway leading from  $\alpha$ -linolenic acid (ALA) to eicosapentaenoic acid (EPA) and docosahexaenoic acid (DHA).<sup>1</sup> Dietary stearidonic acid increases tissue levels of EPA better than does dietary ALA.<sup>2</sup> Stearidonoyl glycine consists of glycine attached at the carboxy terminus.

#### References

- 1. Guil-Guerrero, J.L. Stearidonic acid (18:4n-3): Metabolism, nutritional importance, medical uses and natural sources. Eur. J. Lipid. Sci. Technol. 109(12), 1226-1236 (2007).
- 2. James, M.J., Ursin, V.M., and Cleland, L.G. Metabolism of stearidonic acid in human subjects: Comparison with the metabolism of other w-3 fatty acids. Am. J. Clin. Nutr. 77(5), 1140-1145 (2003).

WARNING THIS PRODUCT IS FOR RESEARCH ONLY - NOT FOR HUMAN OR VETERINARY DIAGNOSTIC OR THERAPEUTIC USE.

#### SAFFTY DATA

This material should be considered hazardous until further information becomes available. Do not ingest, inhale, get in eyes, on skin, or on clothing. Wash thoroughly after handling. Before use, the user must review the complete Safety Data Sheet, which has been sent via email to your institution.

#### WARRANTY AND LIMITATION OF REMEDY

uyer agrees to purchase the material subject to Cayman's Terms and Conditions. Complete Terms and Conditions including Warranty and Limitation of Liability information can be found on our website.

Copyright Cayman Chemical Company, 06/17/2024

#### CAYMAN CHEMICAL

1180 EAST ELLSWORTH RD ANN ARBOR, MI 48108 · USA PHONE: [800] 364-9897 [734] 971-3335 FAX: [734] 971-3640 CUSTSERV@CAYMANCHEM.COM WWW.CAYMANCHEM.COM