

## Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



## Lieferung & Zahlungsart

siehe unsere Liefer- und Versandbedingungen

## Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

### SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

linkedin.com/company/szaboscandic in



# Agrisera

This product is **for research use only** (not for diagnostic or therapeutic use)

contact: support@agrisera.com

Agrisera AB | Box 57 | SE-91121 Vännäs | Sweden | +46 (0)935 33 000 | www.agrisera.com

#### Product no AS18 MacroStab-100

#### Agrisera Universal Stabilizer (100 ml)

#### **Product information**

Agrisera HRP and ALP-conjugate and protein stabilizer preserves and stabilizes the activity of HRP or ALP conjugates and proteins and other bio-molecules used in following assays: western blot, ELISA, Southern/Northern Blots and lateral flow strips. It is designed as a ready to use stabilizer, buffer, diluent and possible coat block for immunoassays.

This stabilizer does not contain mammalian proteins like bovine proteins of BSA. This allows to mnimize background issues, unwanted protein interactions and concerns with possible transmission of bovine spongiform encephalopathy (BSE). Stabilizer is not affected by freezing and is supplied as one component liquid, ready to use formula.

#### **Features**


- no mammalian proteins included
- minimizes background interactions
- stable through freeze cycles (does not contain glycerol)
- can be used to coat plates and as a blocking agent, potentially enabling both blocking and stabilization at the same time
- can be used directly as an assay buffer for high concentration stock and dilute conjugate solutions
- will stabilize diluted conjugates for months at room temperature
- produced based on SureGuard chemical technology
- can be used directly as an assay buffer
- glycerol will not adversely affect the performance of this stabilizer
- antimicrobial agents are inluded to inhibit bacterial growth

Format Liquid

Storage

In the dark at 2°C to 8°C. This product is stable for a minimum of 2 years at 2°C to 8°C. Room temperature applied for extended periods during use and freezing will not adversely affect performance of this product.

#### **Application information**



Comparison between Agrisera HRP-conjugate stabilizer and competitor