

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere Liefer- und Versandbedingungen

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

linkedin.com/company/szaboscandic in





DATA SHEET

RESEARCH USE ONLY www.cellgs.com

PPH314 PODS® Human Gastrin

Description

The product contains the polyhedrin protein co-crystalized with Human Gastrin. Gastrin is an endogenous peptide produced in the stomach that acts as a selective CCK2 receptor agonist and stimulates gastric acid secretion.

Length 62 aa

Molecular Weight 7.3 kDa

Source Spodoptera frugiperda (Sf9) cell culture

Accession Number

Usage Recommendation

PODS® co-crystals provide a depot of proteins which are steadily secreted. It has been estimated that the biological activity of 50 million PODS® co-crystals generates the same peak dose as 3.3 µg of standard recombinant protein. However, at 5 days following the start of seeding the PODS® co-crystals, there are more than 50% of these peak levels still present in the culture system. Ultimately, the amount of PODS® co-crystals that is optimal for a particular experiment should be determined empirically. Based on previous data, we suggest using 50 million PODS® co-crystals in place of 3.3 µg of standard growth factor as a starting point."To control for cross-reactivity with cells or as a negative control, we recommend using PODS® growth factors alongside PODS® Empty crystals, as the latter do not contain or release cargo protein.

Specifications

Alternative Names CCK2 receptor agonist

Endotoxin Level <0.06 EU/ml as measured by gel clot LAL assay

Formulation PODS® were lyophilized from a volatile solution

AA Sequence MADVAGTSNR DFRGREQRLF NSEQYNYNNS KNSRPSTSLY KKAGFQGPWL EEEEEAYGWM DF

Preparation and Storage

Reconstitution PODS® co-crystals may be reconstituted at 200 million co-crystals/ml in water. 20% glucose has a

buoyant density closer to PODS® co-crystals and can be useful for aliquoting.PODS® co-crystals are

highly stable when stored in aqueous solution (pH range 6 - 8).

Stability and Storage Upon receipt, store at 4°C. PODS® co-crystals are stable for at least 1 year when dry and 6 months

when resuspended.