

Produktinformation



Forschungsprodukte & Biochemikalien



Zellkultur & Verbrauchsmaterial



Diagnostik & molekulare Diagnostik



Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



Lieferung & Zahlungsart

siehe unsere Liefer- und Versandbedingungen

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien

T. +43(0)1 489 3961-0

F. +43(0)1 489 3961-7

mail@szabo-scandic.com

www.szabo-scandic.com

linkedin.com/company/szaboscandic in





Product Data Sheet

Tivozanib

Target:

Cat. No.:HY-10977CAS No.:475108-18-0Molecular Formula: $C_{22}H_{19}ClN_4O_5$ Molecular Weight:454.86

Pathway: Protein Tyrosine Kinase/RTK

VEGFR

Storage: Powder -20°C 3 years

4°C 2 years In solvent -80°C 2 years

-20°C 1 year

SOLVENT & SOLUBILITY

In Vitro

DMSO: 25 mg/mL (54.96 mM; Need ultrasonic)

Preparing Stock Solutions	Solvent Mass Concentration	1 mg	5 mg	10 mg
	1 mM	2.1985 mL	10.9924 mL	21.9848 mL
	5 mM	0.4397 mL	2.1985 mL	4.3970 mL
	10 mM	0.2198 mL	1.0992 mL	2.1985 mL

Please refer to the solubility information to select the appropriate solvent.

In Vivo

- 1. Add each solvent one by one: 10% DMSO >> 40% PEG300 >> 5% Tween-80 >> 45% saline Solubility: ≥ 2.5 mg/mL (5.50 mM); Clear solution
- 2. Add each solvent one by one: 10% DMSO >> 90% (20% SBE-β-CD in saline) Solubility: ≥ 2.5 mg/mL (5.50 mM); Clear solution

BIOLOGICAL ACTIVITY

Tivozanib (AV-951; KRN951) is a selective and orally active VEGFR tyrosine kinase inhibitor with IC₅₀ of 0.21, 0.16, 0.24 nM for VEGFR-1, VEGFR-2, VEGFR-3, respectively. Tivozanib inhibits angiogenesis and vascular permeability in tumor tissues and shows antitumor activity. Tivozanib has the potential for the research of metastatic renal cell carcinoma (RCC) [1][2][3].

In Vitro Tivozanib inhibits the phosphorylation of VEGFR-1, VEGFR-2 and VEGFR-3^[2].

MCE has not independently confirmed the accuracy of these methods. They are for reference only.

In Vivo

Tivozanib (1 mg/kg; p.o.; 14 days) suppresses the development of CNV lesions and leds to a significant regression of established CNV, reducing the affected areas by $67.7\%^{[4]}$.

MCE has not independently confirmed the accuracy of these methods. They are for reference only.

PROTOCOL

Kinase Assay

Cell-free kinase assays are done in quadruplicate with 1 μ M ATP to determine the IC₅₀ values of KRN951 against a variety of recombinant receptor and nonreceptor tyrosine kinases^[1].

MCE has not independently confirmed the accuracy of these methods. They are for reference only.

Cell Assay [1]

Cell-based assays are done to determine the ability of KRN951 to inhibit ligand-dependent phosphorylation of receptor tyrosine kinases. Briefly, the cells are starved overnight in appropriate basic medium containing 0.5% fetal bovine serum (FBS). Following the addition of KRN951 or 0.1% DMSO, the cells are incubated for 1 hour and then stimulated with the cognate ligand at 37°C. Receptor phosphorylation is induced for 5 minutes except for VEGFR3 (10 minutes), c-Met (10 minutes), and c-Kit (15 minutes). All the ligands used in the assays are human recombinant proteins, except for VEGF-C, a rat recombinant protein. Following cell lysis, receptors are immunoprecipitated with appropriate antibodies and subjected to immunoblotting with phosphotyrosine. Quantification of the blots and calculation of IC₅₀ values are carried out^[1].

MCE has not independently confirmed the accuracy of these methods. They are for reference only.

Animal Administration [1]

Mice: Cancer cells are s.c. inoculated into the right flank of the athymic rats. Once established, tumors of 1,500 mm³ are surgically excised and smaller tumor fragments (20-30 mg) are s.c. implanted in the right flank of irradiated rats. Oral administration of KRN951 (0.2 or 1 mg/kg) or vehicle is initiated at the day of randomization (day 0). Tumor volume is measured twice weekly with Vernier calipers, and calculated^[1].

MCE has not independently confirmed the accuracy of these methods. They are for reference only.

CUSTOMER VALIDATION

- Sci Transl Med. 2018 Jul 18;10(450):eaaq1093.
- Cancer Cell Int. 2021 Jun 5;21(1):291.
- Pharmaceuticals. 2023, 16(2), 295.
- Technical University of Munich. 24.01.2018.
- Patent. US20170349880A1.

See more customer validations on www.MedChemExpress.com

REFERENCES

- [1]. Motzer RJ, et al. Tivozanib versus sorafenib as initial targeted therapy for patients with metastatic renal cell carcinoma: results from a phase III trial. J Clin Oncol. 2013 Oct 20;31(30):3791-9.
- [2]. De Luca A, et al. Tivozanib, a pan-VEGFR tyrosine kinase inhibitor for the potential treatment of solid tumors. IDrugs. 2010 Sep;13(9):636-45.
- [3]. Eskens FA, det al. Biologic and clinical activity of tivozanib (AV-951, KRN-951), a selective inhibitor of VEGF receptor-1, -2, and -3 tyrosine kinases, in a 4-week-on, 2-week-off schedule in patients with advanced solid tumors. Clin Cancer Res. 2011 Nov 15;17(22):7156-63.
- [4]. Kang S, et al. Antiangiogenic effects of tivozanib, an oral VEGF receptor tyrosine kinase inhibitor, on experimental choroidal neovascularization in mice. Exp Eye Res. 2013 Jul;112:125-33.

 $\label{lem:caution:Product} \textbf{Caution: Product has not been fully validated for medical applications. For research use only.}$

Tel: 609-228-6898 Fax: 609-228-5909

E-mail: tech@MedChemExpress.com

Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA

Page 3 of 3 www.MedChemExpress.com