

Produktinformation



Forschungsprodukte & Biochemikalien
Zellkultur & Verbrauchsmaterial
Diagnostik & molekulare Diagnostik
Laborgeräte & Service

Weitere Information auf den folgenden Seiten! See the following pages for more information!



Lieferung & Zahlungsart siehe unsere Liefer- und Versandbedingungen

Zuschläge

- Mindermengenzuschlag
- Trockeneiszuschlag
- Gefahrgutzuschlag
- Expressversand

SZABO-SCANDIC HandelsgmbH

Quellenstraße 110, A-1100 Wien T. +43(0)1 489 3961-0 F. +43(0)1 489 3961-7 <u>mail@szabo-scandic.com</u> www.szabo-scandic.com

Heparin Lithium salt

MedChemExpress

Cat. No.:	HY-17567B	HQ. O
CAS No.:	9045-22-1	O ^S
Molecular Formula:	(C ₁₄ H ₂₅ NO ₂₀ S ₃)n.xLi	0 10 0 n
Storage:	4°C, sealed storage, away from moisture and light * In solvent : -80°C, 6 months; -20°C, 1 month (sealed storage, away from moisture and light)	HO OH HN S-OH O $O'O$ $O'O$ O XLi

SOLVENT & SOLUBILITY		
In Vitro	H ₂ O : 100 mg/mL (Need ultrasonic)	
In Vivo	1. Add each solvent one by one: PBS Solubility: 100 mg/mL (Infinity mM); Clear solution; Need ultrasonic	

BIOLOGICAL ACTIVITY		
Description	Heparin Lithium salt is an anticoagulant which binds reversibly to antithrombin III (ATIII). Heparin Lithium salt significantly inhibits exosome-cell interactions.	
IC ₅₀ & Target	Antithrombin III ^[1]	
In Vitro	Heparin is a potent anticoagulant drug based on its ability to accelerate the rate at which antithrombin inhibits serine proteases in the blood coagulation cascade. Heparin interactes most tightly with peptides containing a complementary binding site of high positive charge density. Heparin resembles DNA as both are highly charged linear polymers that behave as polyelectrolytes. Heparin is believed to function as an anticoagulant primarily through its interaction with AT III by enhancing AT-III-mediated inhibition of blood coagulation factors, including thrombin and factor Xa. Heparin binds to AT III and thrombin in a ternary complex, increasing the bimolecular rate constant for the inhibition of thrombin by a factor of 2000. Heparin is principally located in the granules of tissue mast cells that are closely associated with the immune response. Heparin makes numerous contacts with both FGF-2 and FGFR-1 stabilizing FGF-FGFR binding. Heparin also makes contacts with the FGFR-1 of the adjacent FGF-FGFR complex, thus seeming to promote FGFR dimerization ^[1] . MCE has not independently confirmed the accuracy of these methods. They are for reference only.	

CUSTOMER VALIDATION

- Nat Commun. 2024 Feb 7;15(1):1150.
- Research Square Preprint. 2023 May 30.

See more customer validations on $\underline{www.MedChemExpress.com}$

Caution: Product has not been fully validated for medical applications. For research use only.

 Tel: 609-228-6898
 Fax: 609-228-5909
 E-mail: tech@MedChemExpress.com

 Address: 1 Deer Park Dr, Suite Q, Monmouth Junction, NJ 08852, USA